

●● dry-sun uv 2.1

Product Description:

The dry-sun uv 2.1 is an electronic drying appliance with additional UV-C light for hygienic preparation. It dries and cares for all types of hearing aids and earmoulds.



Perfect drying by convective flow of warm air



Hygienic preparation by UV-C light for your health



Easy handling by simple touch on a sensor button



Mobility by USB-C standard



Automatic room- and temperature modulation



Automatic turn-off after comfort drying



All features controlled by microprocessor



Highest product quality with 2 year warranty



Energy-efficient, gentle and absolutely safe

Technical data:

- **Drying:** Drying temperature 50 °C (122 °F) +/- 3 °C (37.4 °F) by thermal convection
- **Sensor button:** 1 x touch for 4 hours comfort drying
Touch for 3 seconds for manual switch-off (outside UV-C mode)
- **Switch-off:** Automatic after 4 hours
- **Hygienic preparation:** 2 x 4 minutes UV-C light with 253.7 nm wavelength
Up to 99.9 % of all common germs and bacteria are effectively eliminated
UV-C intensity max. 400 µW/cm²
- **Operation:** 5 V USB-C connector via supplied USB switching power supply or USB 2.0 interface
- **Safety:** Input voltage monitoring at < 4 V and > 6 V
Automatic safety shutdown at temperatures > 70 °C (158 °F)
- **Power consumption:** Average power 2.5 W
Maximum power (UV-C) 4 W
Stand-by power < 0.1 W

Energy supply:

Electronic USB switching power supply (class VI)

- **Input:** 100-240 V, AC 50/60 Hz
- **Output:** 5 V, 1A, DC
- **Type:** EU / UK / US / JP / AUS available

Packaging content:

dry-sun uv 2.1
USB power supply
USB-A to USB-C cable 1.0 m
Hardcover
Instruction for use

Dimensions/Weight:

110 x 90 x 35 mm (L x W x H)
168 g


Warranty:

LOT 2 years

Identification:

WEEE-Reg. no.:
045445-002-1720-1720
UDI-DI 4260277010478

Manufacturer:

 flow-med GmbH
www.flow-med.com



Distribution

audifon GmbH & Co. KG · Werner-von-Siemens-Straße 2 · D-99625 Köllda · Germany
Tel. +49-3635-4056-590 · Fax +49-3635-4056-589 · contact@audifon.com

