

●● multistreamer pro

Bluetooth technology for convenient streaming of external audio sources

- Crystal clear audio from TV
- Enjoy music & podcast from media devices
- Effortless understanding of distant conversation, e.g. at conferences
- Hands-free telephone calls



The audifon multistreamer pro is used to connect your audifon Wings hearing aids to Android and other Bluetooth-enabled devices. It can also be used to transmit an audio signal (e.g. TV signal) or as a table microphone in acoustically demanding situations (e.g. restaurants/ conferences).

multistreamer pro is an optional accessory for all Wings hearing aids from audifon:

- lewi R Li/ risa R Li
- lewi R/ risa R
- lewi S/ risa S

Technical information

- Volume control while music streaming
- Audio streaming (stereo) via A2DP-Control
- Phone connection featuring HFP (handsfree profile) or HSP (headset profile) protocol
- Streaming from 3.5 mm external source audio input (TV/ mp3-player)
- Live-microphone mode
- Lithium-Ion battery with 200 mAh
- Audio bandwidth 7 kHz
- Frequency range: 2402 - 2480 MHz
- Pairing with one or two hearing aids
- Bluetooth 5.0 protocol: Bluetooth Classic and Bluetooth Low Energy

The transferred signals are encrypted and exchanged between the coupled devices.

Rechargeable battery	Size and weight
<ul style="list-style-type: none"> • Standby up to 18 hours • Streaming up to 8 hours • Charging 5 hours 	<ul style="list-style-type: none"> • L 5,8 x W 2,8 x H 1,7 cm • Weight: 17 g
Standards	
EU: <ul style="list-style-type: none"> • EN 300 328 V2.2.2 • EN 301 489-1 V2.2.3 • EN 301 489-17 V3.3.1 • EN IEC 62368-1:2024+A11:2024 	USA: <ul style="list-style-type: none"> • FCC CFR47 part 15.247

CE  FCC number YU2-MS2