

BATTERIES

•• superior quality, longer lasting, environmentaly friendly





audifon's mercury free zinc air batteries offer environmentally friendly and long lasting energy supply for modern hearing systems.

battery size	colour code	capacity	weight
675	blue	650 mAh	approx. 1,83 g
13	orange	300 mAh	approx. 0,83 g
312	brown	170 mAh	approx. 0,58 g
10	yellow	100 mAh	approx. 0,30 g

Technological information

- Available in sizes 675, 13, 312 and 10.
- **audifon** batteries are 100% free of mercury and therefore an environmentally friendly alternative to conventional hearing aid batteries.
- Optimal storage conditions add to the high life period of **audifon** hearing system batteries. Please store at room ambient temperature. Heat and cold may affect life period of batteries.
- The life period of **audifon** mercury free batteries depends on the gain of the respective hearing system and total operation time. **audifon** mercury free zinc air batteries offer a long life period compared to other battery types.



Zentrale

audifon GmbH & Co. KG Werner-von-Siemens-Straße 2 99625 Kölleda Deutschland

Telefon +49-3635-4056-590 Fax +49-3635-4056-589

contact@audifon.com

www.audifon.com