

Solutions for better hearing

u-sonic

# audifon u-sonic





## audifon Solutions for better hearing

Dirt and ear wax can impair the function of ear molds, hearing protection, slim tubes and cause failures. The audifon u-sonic offers optimal cleansing. Powerful ultrasonic at 46 kHz ensures already with clean water hygienic cleanliness. Other applications are small dentures, braces, jewelry, coins and other ultrasonic cleanable parts.

### Product Specification

- all functions are fully processor controlled
- effective cleansing by power-ultrasonic
- ultrasonic frequency at about 46 kHz
- touch sensor button starts 6 min comfort-cleansing program
- hold sensor button >3 seconds for manually stand-by
- automatically stand-by after 6 min cleansing
- tub capacity 100 ml
- short time operation KB10
- protection class IP43
- appliance class II

### Power

- power supply 220–240 V, AC 50–60 Hz
- EU power plug with 1,5 m cord
- power consumption  
max. power: 25 W  
stand-by: <250 mW

### Size

- L 100 mm x W 100 mm x H 84 mm

### Weight

- net: 450 g

### Scope of Delivery

- u-sonic, power plug
- cleansing bag
- user guide

### Warranty

2 years after delivery

### Distribution

audifon GmbH & Co. KG  
Werner-von-Siemens-Straße 2  
99625 Kölleda  
Germany

Phone +49-3635-4056-590  
Fax +49-3635-4056-589

contact@audifon.com



Manufacturer:  
flow-med GmbH · Silberstraße 4 · 71069 Sindelfingen · Germany  
Phone +49-7031-6651-187 · info@flow-med.com · www.flow-med.de