

Solutions for better hearing

dry go

audifon dry go





audifon Solutions for better hearing

audifon's dry go is an easy to use drying station for all kinds of electronic devices. dry go has a USB power supply and is therefore the perfect choice when travelling.

Product Specification

- all functions are fully processor controlled
- effective drying by convective flow of warm and dry air
- two-circute heating with automatic room- and temperature modulation
- drying temperature range about $50^{\circ}\text{C} \pm 3^{\circ}\text{C}$
- touch sensor button: 4 hour comfort-drying starts
- hold sensor button for more than 3 seconds: stop program manually
- programm finishes automatically after 4 hours
- safety interruption to avoid temperatures beyond 70°C
- operation based on 5V USB via supplied USB power supply or USB 2.0 port
- low energy consumption: max. power 2.5W, stand-by: $<0.1\text{W}$

Power

- Switching power supply (class V) with USB-socket and USB cable (1 m)
- Input: 100–240V, AC 50/60 Hz
- Output: 5V, DC, 1A, Adapter EU/UK/US/AUS available upon request

Size

- H 37 mm x W 76 mm x D 92 mm

Weight

- net: 147g

Scope of Delivery

- dry go, USB power supply, USB cable
- travel bag
- user guide

Warranty

2 years after delivery

Distribution

audifon GmbH & Co. KG
Werner-von-Siemens-Straße 2
99625 Kölleda
Germany

Phone +49-3635-4056-590

Fax +49-3635-4056-589

contact@audifon.com



Manufacturer:

flow-med GmbH · Neckarweg 7 · 71139 Ehningen · Germany

Phone +49 7034-2774-718 · info@flow-med.com · www.flow-med.de